

Disclaimer: As MacPractice is not the manufacturer or developer of 3rd party products or services, MacPractice recommends that clients independently verify 3rd party current and ongoing macOS and MacPractice compatibility of 3rd party products with the manufacturer or developer.

	MacPractice Gen 11		MacPractice Gen 12	
	High Sierra (10.13)	Mojave (10.14)	Mojave (10.14)	Catalina (10.15)
3rd Party Applications				
Microsoft Office 2011	Yes	Yes	Yes	No
Microsoft Office 2016	Yes	Yes	Yes	Yes
Microsoft Office 2019	Yes	Yes	Yes	Yes
Microsoft Office 365	Yes	Yes	Yes	Yes
QuickBooks	Yes (QuickBooks 2016+)	Yes (QuickBooks 2016+)	Yes (QuickBooks 2016+)	Yes (QuickBooks 2019+)
Quicken	Yes (Quicken 2016+)	Yes (Quicken 2016+)	Yes (Quicken 2016+)	Yes (Quicken 2016+)
Parallels	Yes (Parallels 13+, Windows XP - 10)	Yes (Parallels 13+, Windows XP - 10)	Yes (Parallels 13+, Windows XP - 10)	Yes, (Parallels 14, Windows XP - 10)
VMWare Fusion	Yes (Fusion 8.5, Windows XP - 10)	Yes (Fusion 11+, Windows XP - 10)	Yes (Fusion 11+, Windows XP - 10)	Yes (Fusion 11+, Windows XP - 10)
MacPractice Services				
AutoRemind	Yes	Yes	Yes	Yes
CDA iTrans	Yes	Yes	Yes	Yes
Change Healthcare	Yes	Yes	Yes	Yes
eStatements (Change Healthcare)	Yes	Yes	Yes	Yes
Integrated Network Faxing	Yes	Yes	Yes	Yes
Online Patient Registration	Yes	Yes	Yes	Yes
Online Patient Scheduling	Yes	Yes	Yes	Yes
Patient Portal	Yes	Yes	Yes	Yes
Reputation Marketing	Yes	Yes	Yes	Yes
Secure Messaging	Yes	Yes	Yes	Yes
Transworld Collections	Yes	Yes	Yes	Yes
TSYS (TransFirst Credit Card Interface)	Yes	Yes	Yes	Yes
AuthPay X	Yes	No	No	No
DollyDrive Internet Backup	Yes	Yes	Yes	Awaiting Vendor Response
X-Ray Native Interface: Native interfaces send images directly to the Digital Radiography Ability in MacPractice				
Native interfaces send images directly to the Digital Radiography Ability in MacPractice				
Acteon SOPIX2 USB Sensor (1, 2)	Yes	Yes	Yes	Yes
Acteon PSPIX Phosphor Plate (0, 1, 2, 3, 4)	Yes	Yes	Yes	Yes
DEXIS Products with DTX Studio	No, this is new in MP Gen 12	No, this is new in MP Gen 12	Yes	Yes
Gendex via DTX Studio	No, this is new in MP Gen 12	No, this is new in MP Gen 12	Yes	Yes
MyRay X-pod with Sensor (1, 2)	Yes (Bluetooth)	Yes (Bluetooth)	Yes (Bluetooth)	Yes (Bluetooth)
Planmeca Prosensor (0, 1, 2)	Yes	Yes	Yes	Yes
Planmeca Romexis	Yes	Yes	Yes	Yes
QuickRay USB Sensor (1, 2)	Yes	Yes	Yes	Yes
SuniRay/SuniRay 2 Sensor (1, 2)	Yes	Yes (Unsupported)	Yes (Unsupported)	Yes (Unsupported)
X-Ray Native Interface - Pan, Ceph, Cone Beam				
Gendex via DTX Studio	No, this is new in MP Gen 12	No, this is new in MP Gen 12	Yes	Yes
iCAT via DTX Studio	No, this is new in MP Gen 12	No, this is new in MP Gen 12	Yes	Yes
Instrumentarium via DTX Studio	No, this is new in MP Gen 12	No, this is new in MP Gen 12	Yes	Yes

InvivoMac (Anatomage) - Most Devices	Yes	Yes	Yes	Yes
MyRay Hyperion	Yes	Yes	Yes	Yes
Owandy via DTX Studio	No, this is new in MP Gen 12	No, this is new in MP Gen 12	Yes	Yes
Planmeca ProOne and ProMax	Yes	Yes	Yes	Yes
X-Ray Bridge (Windows)				
Please verify your computer or VM Software can run the appropriate version of Windows that your Digital Radiography Software Requires.				
Apteryx/CliniView (Most Sensors, Pan, Ceph)	Yes	Yes	Yes	Yes
CADI	Yes	Yes	Yes	Yes
DEXIS	Yes	Yes	Yes	Yes
Lightyear	Yes	Yes	Yes	Yes
OpTime (Phosphor Plate)	Yes	Yes	Yes	Yes
ProfSuni	Yes	Yes	Yes	Yes
ScanX / Visix (Phosphor Plate)	Yes	Yes	Yes	Yes
Schick (Pre-Sirona)	Yes	Yes	Yes	Yes
Sirona Sidexis 4	Yes	Yes	Yes	Yes
XDR	Yes	Yes	Yes	Yes
Imaging Software Bridge (Mac)				
Anatomage - InvivoMac	Yes	Yes	Yes	In Development
DTX Studio	No, this is new in Gen 12	No, this is new in Gen 12	Yes	Yes
Planmeca Romexis	Yes	Yes	Yes	Yes
RadioVision	Yes	Yes	Yes	Yes
Perio Charting Data Entry - Foot Pedal				
Dental R.A.T.	Yes	Yes	Yes	Yes
Plug and Play Cameras				
DentSmart Dental Camera	Yes	Yes	Yes	Yes
iCam	Yes	Yes	Yes	Yes
ImageMaster Billet USB Intraoral Camera	Yes	Yes	Yes	Yes
ImageMaster USB	Yes	Yes	Yes	Yes
MouthWatch	Yes	Yes	Yes	Yes
RCI RealCloud HD1 Mac	Yes ‡	Yes ‡	Yes ‡	Yes ‡
‡ Please note, this camera is functional and images can be taken in MacPractice, but the software for the button on the camera is still in development.				
Camera (With Interface Software)*				
These cameras can be used with MacPractice by using interface software (listed below) and saving images to an import folder				
Digital Doc Iris Mac Version	Yes	Yes	Yes	No*
Discovery Cam	Yes	Yes	Yes	No*
Discovery Cam 2	Yes	Yes	Yes	No*
DocPortMacro USB	Yes	Yes	Yes	No*
DR Cam	Yes	Yes	Yes	No*
SoproCare	Yes	Yes	Yes	Yes, with Sopro Imaging MAV2M
SoproLife	Yes	Yes	Yes	Yes, with Sopro Imaging MAV2M
Sopro 617	Yes	Yes	Yes	Yes, with Sopro Imaging MAV2M
Sopro 717	Yes	Yes	Yes	Yes, with Sopro Imaging MAV2M

SUNI SuniCam II USB	Yes	Yes	Yes	No*
Camera Interface Software				
DTX Studio	No, this is new in Gen 12	No, this is new in Gen 12	Yes	Yes
ProView USB Software	Yes	Yes	Yes	No*
SoproCapture	Yes	Yes	Yes	No*
Sopro Imaging MAV2M	Yes	Yes	Yes	Yes
VideoGlide	Yes	Yes	Yes	No*
* The Camera Interface software is not compatible with macOS 10.15 Catalina and higher. Please contact MacPractice if you have any questions about compatibility of Intra Oral cameras.				
Fujitsu ScanSnap				
Fujitsu announced that it would not be releasing Catalina drivers for some models of their ScanSnap scanners. Devices with Catalina drivers: iX1500, iX500, iX100, S1300i, S1100i, S1100, SV600.				
If your scanner is not listed, VueScan allows the unsupported ScanSnap devices to work in Catalina.				
Printers				
Please verify that your laser printer is compatible with the version of macOS that you will be printing from. Most laser printers with postscript emulation will print insurance and other forms from MacPractice, and all printers should print word processing documents.				